KKD[®] provisional CrownResin

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Provisional CrownResin information:

self-curing temporary crown and bridge resin

KKD® provisional CrownResin is a self curing resin for constructing temporary crowns, bridges and splints. It offers many advantages, from the simple technic to first-class results.

Good compatibility with vital stumps and tissue. KKD® provisional CrownResin does not contain methyl methacrylate. Due to the special formula of powder and liquid a perfect polymerization is obtained. There is practically no residual monomer which could effect vital stumps.

Great stability on multi-unit bridges and thin-ending crown rims due to increased hardness of the material. KKD® provisional CrownResin makes it possible to construct even multi-unit bridges on restorations with large stumps and little space.

Easy handling in three steps:

1. Mixing phase: Powder and liquid are mixed to a doughy, guaranteed bubble-free consistency. This avoids air inclusion in the temporary restoration.
2. Plastic phase: Period 3 to 5 minutes. There are no distortions while removing from the mouth due to the extremely high withdrawal of the material.
3. Curing phase: The final polymerization takes place during the next 5 minutes. For faster curing a pressure pot can be used. This also results into an even better color stability and optimum, practically monomerfree polymerization.

Security instructions: The liquid is irritable. It has to be stored at a cool place. Immediately close bottle after every use.

Individuals:

ea. 100 ml powder	universal (20 + 21)*	REF 23214	50 ml liquid	REF 23216
	dark (31 + 32)*	REF 23212	150 ml liquid	REF 23215
	light (12 + 13)*	REF 23213		

*colors correspond to Biodent colors

Presentation:
100 ml powder (color universal),
50 ml liquid, 1 mixing bowl,
1 mixing spatula
REF 232065

Directions:

- 1 Prepare the silicone or alginate impression as usual.
2. To mix provisional CrownResin, pour liquid in a mixing bowl and add sufficient powder in the desired color until a doughy consistency is obtained. This mixture can now be placed immediately into the impression.
3. Insert the impression furnished with KKD® provisional CrownResin into the mouth of the patient and leave it for 3 minutes.
4. KKD® provisional CrownResin has now reached its plastic phase and can be easily removed with the impression from the mouth, even if there are undercuts on the preparation. The withdrawal of the material is excellent.
5. Remove the temporary restoration from the impression and place it for security and to establish a proper occlusion over the prepared stumps. After this control remove it from the mouth again. In this plastic phase (period 3 to 5 minutes) KKD® provisional CrownResin can be easily trimmed with a scissor which simplifies the finishing of the temporary restoration.
6. During the next 5 minutes the final polymerization takes place. For faster curing we recommend to use a pressure pot filled with warm water of approx. 55° C and cure at a pressure of 2 bar. This also results into an even better color stability, optimum (practically monomerfree) polymerization and higher mechanical strength. We recommend TOPDENT Polypot, the practical and safe polymerization unit (independent of compressed air, no electrical parts).
7. Finish the temporary restoration as usual and insert with temporary cement. Temporary restorations with the KKD® provisional CrownResin can of course be constructed after any common technics (e. g. indirect method).

Provisional CrownResin information:

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KENTZLER-KASCHNER DENTAL GmbH