astron®

Exclusive in Europe
KENTZLER-KASCHNER DENTAL GmbH
Mühlgraben 36
73479 Ellwangen/Jagst
Telefon: +49 (7961) 9073-0
Fax: +49 (7961) 52031
e-Mail: info@kkd-topdent.de

astron[®]

Easier for dentists, better for patients

Technology for dental perfection -
inspired by nature

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REF 26178

for comfort nightguards

- extremely break-resistant
- thermoactive flexibility
- free of MMA
- self-adjusting
- always transparent, no yellowing
- easy to polish (like acrylic resins)
- fast and easy processing
- the thermoactive flexibility provides precise adaption
- better adaption due to tension-free wearing comfort
- no feeling of pressure in the mouth

Ceramic mixing bowl

When kept in the fridge for a short time (5-10 min.), it provides a longer processing time of the thermoactive special resin. No release of chemical components, no chemical bond between bowl and CLEARsplint®!

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
**“CLEARsplint is perfect for a premium splint.
Doctors love the material and keep coming back for more.”** Lee, Glidewell Laboratories

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

Instructions

Store powder and liquid in the refrigerator at 4 - 8 °C. (Practical tips: prior to use put mixing bowl into the refrigerator for 30 minutes to extend processing time - we recommend the use of ceramic - or glass bowls to avoid any chemical bond with the material). Thus working time will be extended of this heat-sensitive material. In order to avoid uncontrollable chemical reactions please use only mixing bowls and -spatulas which do not get in contact with other resins. Processing repairs as usual. Prior to processing repairs, put some CLEARsplint® liquid on area to be repaired and let it dry for some seconds.

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1. prepare and measure model and mark tooth equator - splint can be made above tooth equator
2. wax up splint with a wax plate
3. prepare protection wall and let it cure. Remove it and prepare on the tubers, resp. on both sides of the molar triangle an overflow canal
4. boil out the model and water it
5. isolate the still warm plaster model, dry it and fix with protection wall
6. remove powder and liquid from the refrigerator just prior to mixing and pouring
7. mixing ratio: 10 g powder to 7 g liquid (**1 g powder to 0.7 g liquid to or a multiple**)
Important: respect mixing ratio thoroughly! For one splint approx. 8 - 10 g powder are required
8. mix it thoroughly for approx. 5 - 10 seconds as quickly as possible and bubble-free (use separate measuring cup and -spatula)
ATTENTION: pour resin immediately into protection wall and cure it as quick as possible at 2 bar and 50 - 53 °C (122 - 127 °F) for 45 minutes
9. remove denture mold from pressure pot and place it under cold running tap water for 15 minutes before removing splint from pour mold
10. for removing splint from the plaster model, warm it up at approx. 50 °C (122 °F) and remove it carefully
11. trim and polish as usual
12. immerse splint into handwarm water prior to placing it on the model or in the mouth (flexibilisation)

Press packing Instructions

1. prepare and measure model and mark tooth equator - splint can be made above tooth equator
2. wax up splint with a wax plate and invest model into flask
3. prepare counter part and after curing heat flask in water bath
4. remove wax, boil out model and counter part, water and isolate them
5. remove powder and liquid from the refrigerator just prior to mixing
6. mixing ratio: 10 g powder to 7 g liquid (**1 g powder to 0.7 g liquid to or a multiple**)
Important: respect mixing ratio thoroughly! For one splint approx. 8 - 10 g powder are required
7. mix thoroughly for approx. 5 - 10 seconds as quickly as possible and bubble-free (use separate measuring cup and spatula)
8. as soon as mix reaches a doughy state (3 - 5 minutes) pack into regular denture mold using normal press-packing procedure
9. immerse in 50 - 53 °C (122 - 127 °F) warm water and cure for 60 minutes.
10. cool for 15 minutes under cold running tap water before removing from mold
11. trim and polish as usual
12. for removing splint from the plaster model, warm it up at approx. 50 °C (122 °F) and remove it carefully
13. trim and polish as usual
14. immerse splint into handwarm water prior to placing it on the model or in the mouth (flexibilisation)

Care of the astron® CLEARsplint® nightguard: We recommend astron® DuraClean® REF 26170S

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- licence-free
- no high equipment costs
- easy to enlarge, repair and reline
- chemical bond with prefabricated tooth
- metal-free
- transparent

also suited for:

- > long-term temporary restorations
- > gingival masks
- > dentures for children

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Instructions

Processing: (fast, easy, effective)

The prefabricated teeth are set into the gap to be substituted and the temporary restoration is waxed up in the usual manner (modelled in wax). Retention elements can be modelled in the form of a clamp or as „ gum shield“. Make sure that clamps are modelled in an appropriate thickness to ensure that they are entirely poured out.

We recommend that a bigger part of the retention elements lies in regions with undercuts (below tooth equator), in order to guarantee a safe holding of the partial denture after completion.

Once modellation is finished prepare silicone protection wall.

Boil out the model, water it and isolate the still warm model. Prefabricated teeth are sandblasted or grinded, additionally furnished with mechanical retentions and moistened with CLEARsplint® liquid. A chemical bond between the teeth and CLEARsplint® will be established. Therefore the use of a connector is not absolutely necessary (bond also between acrylic resin and CLEARsplint®).

Reset prefabricated teeth back in the protection wall. Using a scalpel make one hole at each end of the protection wall.

Hint:

Put a small funnel, used in the model casting technique into the filling hole to create a small reservoir, out of which the resin may flow.

Once protection wall has been put on the model, place it for approx. 30 minutes in the refrigerator for cooling. (User tips: prior to use put mixing bowl into the refrigerator for 30 minutes to extend processing time - we recommend the use of ceramic- or glass bowls to avoid any chemical bond with the material). Mix the appropriate (see instructions of use) cool liquid- and powder amount of CLEARsplint® and pour resin immediately into protection wall.

Cure model for 45 minutes at 50 - 53 °C in the pressure pot. Cool model acc. to instructions of use.

For removing restoration from the plaster model, warm model up at approx. 50 °C and remove it carefully. Trim and polish as usual. Your result will be a thermoactive, MMA-free, reliable and enlargable partial denture.

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The premium nightguard material especially for the digital manufacturing of functional therapeutic „comfort splints“ with the CAD/CAM technique. Clearsplint® CAD/CAM Pro possesses the same outstanding advantages of CLEARsplint® as e.g. the thermoactive flexibility and the extreme break resistance.

Some individual work steps can be saved thus reducing the working time enormously. Please use Clearsplint® CAD/CAM Pro only in blanks out of the Starter Kit or in device-specific appropriated transparent blanks.

The results are marked by closest fit and a surface which has to be simply high gloss polished.

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Ø 98x20 mm with shoulder
For all big milling units
e.g. Wissner, Röders, Datron, imes-icore, vhf, Roland as well as Sirona inLab MC X5

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Instructions

Store powder and monomer in the refrigerator at 4 - 8 °C. (Practical tips: prior to use put mixing bowl into the refrigerator for 30 minutes to extend processing time - we recommend the use of ceramic - or glass bowls to avoid any chemical bond with the material). Thus working time will be extended of this heat-sensitive material. In order to avoid uncontrollable chemical reactions please use only mixing bowls and -spatulas which do not get in contact with other resins. Processing repairs as usual. Prior to processing repairs, put some Clearsplint® CAD/CAM Pro liquid on area to be repaired and let it dry for some seconds.

ATTENTION:

Please respect mixing ratio, polymerisation temperature, -pressure and -time thoroughly!
Otherwise material risks not to cure properly or product-specific advantages may be influenced negatively!
• choice of the appropriate position and size of the area to be milled necessary for the filling with Clearsplint® CAD/CAM Pro in the device-specific blank (base blank)

ATTENTION:

Only use blanks which are included in the Starter Kit or device-specific appropriate transparent blanks recommended for the manufacture of splints! Please make sure that the Clearsplint® CAD/CAM Pro resin is well bonded to the base blank! A blank can be refilled repeatedly since material adheres to the residual material of the last filling (Reparability).

Mixing ratio: 1 g powder to 0.7 g monomer

CREATE THE BLANK:

For one splint you need according to size of the milled-out space approx. 30 to 40 g powder. This means you need for **30 g (0.07 lb) powder 21 g (0.05 lb) monomer, resp. for 40 g (0.09 lb) powder 28 g (0.06 lb) monomer.**


- Isolate base plate with a thin layer of vaseline
- fix milled blank with the clamps onto the Clearsplint® CAD/CAM Pro base (base plate)
- take powder and monomer out of the refrigerator just immediately prior to use. Mix thoroughly the appropriate amount for approx. 5 - 10 seconds as quick as possible and bubble-free
- fill the milled-out space as quickly as possible with the Clearsplint® CAD/CAM Pro resin and wait until resin forms a membrane (start of the surface polymerisation)
- This is absolutely necessary to avoid water inclusions during polymerisation!
- Immerse the filled and fixed blank in the polymerisation pot. Make sure that the complete blank is covered with water
- polymerisation in 50 - 53 °C (122 - 127 °F) warm water in a pressure pot at 2 bar (29 PSI) for 45 minutes
- remove Clearsplint® CAD/CAM Pro base from pressure pot and let it cool down under cold running tap water
- then clamp blank into the milling unit and start milling process of your Clearsplint® CAD/CAM Pro splints made of Clearsplint® CAD/CAM Pro resin can be trimmed and polished as usual

HINT:

We recommend to water the finished splint for approx. 3 - 4 hours, then to rinse and brush it under running water with a soft brush and to disinfect it afterwards with an alcohol-free disinfectant.

ATTENTION:

Please write down the lot number of the used resin, as well as the date of the manufacturing and the name of the patient. For the care of the splint we recommend to use DuraClean², a cleaner for splints especially developed from the astron company.
Do not use denture cleaners, as well as abrasive tooth pastes, since they may damage the splints even up to complete destruction of the Clearsplint® CAD/CAM Pro „comfort splint“.

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
- industrially manufactured disc made of CLEARsplint®
- no mixing and filling anymore
- constantly high quality standard (almost bubble-free)
- enormous time-savings, use without time-consuming prep work
- **MMA-free**

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- easy to handle
- no need for special devices
- **licence-free**
- natural colour
- excellent wearing comfort
- very break-resistant
- easy to relin
- suitable for all manufactured teeth

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
“Astron is the best material I have ever used period...nothing even comes close and I’ve tried everything.” Dan, Flexease D/L

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Instructions - Labtec Pro

Store powder and liquid in the refrigerator at 4 - 8 °C. (Practical tips: prior to use put mixing bowl into the refrigerator for 30 minutes to extend processing time – we recommend the use of ceramic – or glass bowls to avoid any chemical bond with the material). In order to avoid uncontrollable chemical reactions please use only mixing bowls and - spatulas which do not get in contact with other resins. For RELININGS and REPAIRS we recommend to use Labtec Reline.

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1. prepare model and wax up denture in regular/usual manner (use natural stone class III)
2. prepare protection wall and let it cure
3. remove it and prepare on the tubers, resp. on both sides of the molar triangle two overflow canals
4. boil plaster model out and water it
5. isolate the still warm plaster model
6. remove the isolation immediately under running tap water and dry it
7. sandblast the teeth, make retentions and reset the teeth back in the protection wall
8. fix the model with the protection wall and put it into the refrigerator for 15 minutes
9. powder and liquid have to be precisely-weighed and then put the appropriate liquid- and powder amount for approx. 15 minutes in the refrigerator
10. **Mixing ratio: ideally, use the measuring spoons included in the package**, or 1.0 g powder to 0.85 g liquid or a multiple, using a cold ceramic bowl (for one denture 20 g powder to 17 g liquid are required), spatulate slowly and bubble-free
11. pour resin immediately in a thin stream into the protection wall, immerse in pressure pot and cure for 60 minutes at a temperature of 52°C and 2 bar
12. remove denture from pressure pot and prior to open the unit, cool it for 15 minutes under cold running tap water
13. after reocclusion remove denture from plaster model and finish and polish the denture as usual

Attention: Do not steam off the denture, otherwise the high-gloss of the surface will be damaged and becomes matt. Ultrasonic cleaning and rinsing with fresh water and detergent are sufficient.

Hint: for easier removal of the denture from the plaster model, put model without „protection wall“ for approx. 30 seconds in warm water and remove it carefully

Press packing Instructions

1. prepare model and wax up denture in regular/usual manner (use natural stone class III)
2. invest model into flask
3. prepare counter part and after curing heat flask in water bath
4. open the flask, remove wax, boil out model and counter part, water and isolate them
5. sandblast the teeth, make retentions and reset the teeth back in the flask
6. put the flasks into the refrigerator for 15 minutes
7. remove powder and liquid from the refrigerator just prior to mixing
8. **Mixing ratio: ideally, use the measuring spoons included in the package**, or 1.0 g powder to 0.85 g liquid or a multiple, using a cold ceramic bowl (for one denture 30 g powder to 25,5 g liquid are required), spatulate slowly and bubble-free
9. important: as soon as mix reaches a doughy state (5 – 8 minutes) pack into regular denture mold using normal press-packing procedure
10. close mold slowly to permit resin to flow. Place clamp and remove it together with flask
11. immerse flask in 52°C warm water and cure for 2 hours
12. cool flask for minimum 15 minutes under cold running tap water before removing denture from mold
13. after reocclusion remove denture from plaster model and finish and polish the denture as usual

Attention: Do not steam off the denture, otherwise the high-gloss of the surface will be damaged and becomes matt. Ultrasonic cleaning and rinsing with fresh water and detergent are sufficient.

Hint: for easier removal of the denture from the plaster model, put model for approx. 30 seconds in warm water and remove it carefully

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Relining Unit »Okklamat«

For a fast relining of partial and complete dentures. Suitable for any polymerisation unit.

REF 8081010

astron Labtec Reline is a resin for relinings, repairs and enlargements in the partial denture- and full denture prosthetics. It is anallergenic, free of MMA and ensures optimal bonding with the base. Universal for use - no special bonding agent required.

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Ceramic mixing bowl

When kept in the fridge for a short time (5-10 min.), it provides a longer processing time of the special resin. No release of chemical components, no chemical bond between bowl and resin!

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
Instructions


Relining

1. Place model together with the denture which has to be relined in the relining unit and isolate the denture as usual for counter mold preparation
 2. After curing of the stone counter mold dismantle the unit and remove carefully denture from the model
 3. Remove impression material from the denture and roughen the denture which has to be relined – e.g. with sandpaper
 4. Boil plaster model out and water it
 5. Isolate the still warm plaster model
 6. Remove the isolation immediately under running tap water and dry it
 7. Reposition denture in the counter mold and fix it with e.g. sticky wax
 8. Brush Labtec liquid as bonding agent onto the roughened denture
 9. Powder and liquid have to be precisely-weighed and then put the appropriate liquid and powder amount for approx. 15 minutes in the refrigerator
 10. **Mixing ratio: ideally, use the measuring spoons included in the package,** or 1.0 g powder to 0.85 g liquid or a multiple, using a cold ceramic bowl
 11. Spatulate powder slowly and bubble-free into liquide
 12. Working time just 2 - 3 minutes
 13. Place counter mold with denture in the relining unit
 14. Immerse in pressure pot and cure for 60 minutes at a temperature of 52 °C and 2 bar
 15. Remove denture mold from pressure pot and prior to open the unit, cool it for 15 minutes under cold running tap water
 16. Now open carefully the relining unit, remove denture from plaster model and finish and polish the denture as usual
- Attention: Do not steam off the denture, otherwise the high-gloss of the surface will be damaged and becomes matt. Ultrasonic cleaning and rinsing with fresh water and detergent are sufficient**
- Hint: for easier removal of the denture from the plaster model, put model for approx. 30 seconds in warm water and remove it carefully**

Enlargement and repair

1. Prepare model in regular manner, prepare denture and, if necessary, set up teeth in wax
 2. Prepare „protection wall“ and boil it out
 3. Roughen denture resin and, if necessary, sandblast the tooth and make retentions with a bur or drill holes
 4. Apply Labtec liquid as bonding on the area to be repaired
 5. **Mixing ratio: ideally, use the measuring spoons included in the package,** or 1.0 g powder to 0.85 g liquid or multiple, using a cold ceramic bowl
 6. Spatulate slowly and bubble-free
 7. Pour resin immediately in a thin stream into the protection wall, immerse in pressure pot and cure for 60 minutes at a temperature of 52 °C and 2 bar
 8. Remove denture mold from pressure pot and cool it for 15 minutes under cold running tap water
 9. Remove denture from protection wall and finish and polish the denture as usual
- Attention: Do not steam off the denture, otherwise the high-gloss of the surface will be damaged and becomes matt. Ultrasonic cleaning and rinsing with fresh water and detergent are sufficient**
- Hint: for easier removal of the denture from the plaster model, put model without „protection wall“ for approx. 30 seconds in warm water and remove it carefully.**

DuraClean²TM

"Finding a decent nightguard cleaner at the drugstore was impossible! Thanks for coming out with DuraClean²TM ... it really works!" Cheryl, Las Vegas, NV

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Special cleaning agent for
CLEARsplint® nightguards and
interims, as well as Labtec denture
materials

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**Cleaning with a toothbrush and
toothpaste**

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after

After a short period of time
the surface is already scratched and matt

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after

Even after repeated cleaning, surface stays clear and
undamaged

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Due to the high abrasive effect of tooth pastes their
permanent use may result in optical impairments (grey
shadows) of your highly aesthetic transparent CLEARsplint®
premium splint.

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
Aggressive cleaning agents, such as denture cleaners must
not be used with CLEARsplint® premium splint at all.

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set cleaning powder (170 g) with box and brush
Eco-kit (6 sets)
cleaning powder
Eco-Kit cleaning powder (6 tins)
cleaning brush
cleaning box

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Special burs for the preparation of relinings and dimensional changings of soft and hard denture resins. Well-balanced, long-lasting; produce a rough surface for optimal retention.

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

Allows an easy cleaning of 14 acrylic teeth by steam jet cleaner or wax scalding unit. With the tooth arrangement mixing teeth up is nearly eliminated. Due to the special colour codings the sieves can be easily assigned to specific works. Thermostable plastic material.

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ca./ approx. 225 g
 REF **26127**

TD Orange Wheel - medium Pink Wheel - soft

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Impregnated polishing wheels for high lustre on acrylic surfaces, especially on thermoplastics and soft relinings.

- High-lustre without melting effects

⊙ 5.000 - 10.000 min⁻¹

Körnung: mittel/ Grit: medium
 22 x 3 mm - 10 St./ pcs. REF **22459**

Körnung: fein/ Grit: fine
 22 x 3 mm - 10 St./ pcs. REF **22458**

astron®

EC REP

Exclusive in Europe

KENTZLER-KASCHNER DENTAL GmbH
Mühlgraben 36
73479 Ellwangen / Jagst
Telefon: +49 (7961) 9073-0
Fax: +49 (7961) 52031
e-Mail: info@kkd-topdent.de