YDM PARODONTOMETER

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
#5 (PCP11) REF 12460

CPITN-C (P-WHO) REF 12461

READ ALL INSTRUCTIONS BEFORE PROCEEDING

CAUTION: Read following instructions for use.

Gram Probe is designed to measure periodontal pocket. Please do not use it for other purposes.

An error may occur when the strength is put in narrow position.

Loading 20 to 25g, the pin hits part A. Please do not put further weight.

Please do not touch movable part B (the black part) during the measuring of periodontal pocket.

Do not turn the screw C.

HANDLING AFTER USE

Please remove foreign substances (chemical, blood, humors and tissues) adhering to the product as soon as possible. It becomes the cause of corrosion. We recommend glutaraldehyde as disinfectants.

Do not use steel tools such as wire brush to clean.

DISINFECTING & STERILIZING

The tip with black indicated measurements may peel off when having contact with other instruments. We recommend sterilizing it separately in autoclavable pouch or wrapping it with a gauze.

When using disinfectants, please read the directions carefully.

STORAGE

Fully dry before you store the product. When drying is inadequate, it may lead to corrosion.

Please do not store the product together with instruments with rust or chemicals.

STERILIZATION AND DISINFECTING

1. AFTER USING THE INSTRUMENT, RINSE WITH RUNNING WATER.

2. ULTRASONIC CLEANING AND WASH WITH DISINFECTANT.

- please remove foreign substances (chemical, blood, humors, tissues, etc.) from the instrument before going through process of sterilization to prevent corrosion and decrease in the effect of sterilization.
- adhesion of blood or protein can effectively be removed by using disinfectants.
- do not place mirror and periodontal probe in ultrasonic cleaning.
- do not use steel tools such as cleanser, wire brush, steel wool and file to remove foreign substances.
- do not use home detergent, it will likely to result in corrosion of metal.

3. RINSE WITH WATER

- rinse instruments with water and dry before sterilization.

4. STERILIZE

- avoid autoclaving that uses chemical gas and solutions such as alcohol
- we recommend glutaraldehyde as cold sterilizing solution. Do not expose YDM products in Sodium Hypochlorite, Benzalkonium Chloride, Benzethonium Chloride, Providone-Iodine, Formalin, Phenol, Chlorhexidinegluconate.
- do not use super acidic water.
- use purification water. Tap water may lead to corrosion.

5. DRY & STORAGE

- after going through the process of cold sterilization, please rinse well with water and well dry before stocking.
- joint of extracting forceps, pliers, scissors, hemostatic forceps and needle holders are recommended to be lubricated periodically with a siliconebased lubricant, p.ex. KKD Silikone Spray. This will avoid rust and maintains smooth movement of an instrument.
- please do not stock instruments with others that are already corroded.
- please do not stock instruments with any chemical substances. It may lead to corrosion.