

# TOPDENT PERLADIA

Schleifgummi mit großem Füllgehalt  
hochreiner Diamanten

NEU

Durch die einzigartige Klebtechnik  
von Gummi und Diamanten erreicht  
der PERLADIA gleichzeitiges

**Abtragen und Polieren von Zirkon,  
Presskeramik und Keramik in einem  
Arbeitsgang** auf Nanoebene.

Die Plaqueeanhftungen an keramischen  
Restorationen werden dadurch sehr  
gering gehalten.

- Lange Lebensdauer
- hochwertige Qualität und Verarbeitung
- Nutzungsbereich bis zum Kern (Mandrell)
- effizienter Abtrag bei gleichzeitiger Nano-Politur
- geringe Wärmeentwicklung durch  
einen hohen Diamantanteil



Linse dunkelbraun  
Ø 25 mm  
REF 20630



Spitze dunkelbraun  
Ø 5 mm, L 12 mm  
REF 20631

Linse rosa  
Ø 25 mm  
REF 20650



Spitze rosa  
Ø 5 mm, L 12 mm  
REF 20651

**dunkelbraun = mittlere Körnung**  
Hoher Abtrag bei gleichzeitiger Politur

**rosa = feine Körnung**  
Geringer Abtrag bei gleichzeitiger Hochglanzpolitur

# TOPDENT PERLADIA

Abrasive rubbers with large filling content of high pure diamonds

NEW

Due to the unique bonding technology of rubber and diamonds the PERLADIA enables simultaneous **removal and polishing of zirconia, pressed ceramics and ceramics in one operation** on nano level.

Therefore plaque adhesions to ceramic restorations are kept to a minimum.

- Long service life
- Can be used up to the head of the mandrel
- High quality and workmanship
- Efficient removal with simultaneous nano polishing
- Low heat generation due to high diamond content



Lens dark brown  
Ø 25 mm  
REF 20630



Point dark brown  
Ø 5 mm, L 12 mm  
REF 20631

Lens pink  
Ø 25 mm  
REF 20650



Point pink  
Ø 5 mm, L 12 mm  
REF 20651

**dark brown = medium grit**

High removal with simultaneous polishing

**pink = fine grit**

Low removal with simultaneous high gloss polishing